

Product datasheet for **HK200143**

ERCC8 Human qPCR Template Standard (BC009793)

Product data:

Product Type:	qPCR Template Standards
ACCN:	BC009793
Synonyms:	CKN1; Cockayne syndrome WD-repeat protein CSA; CSA; DNA excision repair protein ERCC-8; excision repair cross-complementing rodent repair deficiency, complementation group 8
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP200068) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)