

Product datasheet for **HK200144**

EWSR1 Human qPCR Template Standard (BC000527)

Product data:

Product Type:	qPCR Template Standards
ACCN:	BC000527
Synonyms:	bK984G1.4 (Ewing sarcoma breakpoint region 1 protein); Ewing sarcoma breakpoint region 1; Ewing sarcoma EWS-Fli1 (type 1) oncogene; EWS
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP208415) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)