

Product datasheet for **HK200145**

FAM125A (MVB12A) Human qPCR Template Standard (BC011840)

Product data:

Product Type:	qPCR Template Standards
ACCN:	BC011840
Synonyms:	CFBP; CIN85/CD2AP family binding protein; family with sequence similarity 125, member A; FLJ32495
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP216960) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)