

Product datasheet for **HK200149**

FAM19A2 Human qPCR Template Standard (BC028403)

Product data:

Product Type:	qPCR Template Standards
ACCN:	BC028403
Synonyms:	DKFZp761E1217; DKFZp781P0552; family with sequence similarity 19 (chemokine (C-C motif)-like), member A2; MGC42403; TAFA-2; TAFA2
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP218865) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)