

Product datasheet for **HK200154**

FCHSD2 Human qPCR Template Standard (BC010394)

Product data:

Product Type:	qPCR Template Standards
ACCN:	BC010394
Synonyms:	FCH and double SH3 domains 2; KIAA0769; nervous wreck homolog; NWK; SH3 multiple domains 3; SH3MD3
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP211179) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)