

Product datasheet for **HK200157**

FIP1L1 Human qPCR Template Standard (BC017724)

Product data:

Product Type:	qPCR Template Standards
ACCN:	BC017724
Synonyms:	DKFZp586K0717; FIP1 like 1; FIP1 like 1 (S. cerevisiae); FLJ33619; rearranged in hypereosinophilia; Rhe
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP215465) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)