

## Product datasheet for **HK200164**

### FUBP1 Human qPCR Template Standard (BC017247)

#### Product data:

Product Type:	qPCR Template Standards
ACCN:	BC017247
Synonyms:	DNA helicase V; far upstream element (FUSE) binding protein 1; far upstream element-binding protein; far upstream element binding protein; FBP; FUBP; FUSE-binding protein; OTTHUMP00000038483
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP207343) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 <sup>7</sup> copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)