

Product datasheet for **HK200169**

GALE Human qPCR Template Standard (BC001273)

Product data:

Product Type:	qPCR Template Standards
ACCN:	BC001273
Synonyms:	FLJ95174; FLJ97302; galactose-4-epimerase, UDP-; galactowaldenase; OTTHUMP00000002991; OTTHUMP00000044857; SDR1E1; short chain dehydrogenase/reductase family 1E, member 1; UDP-galactose-4-epimerase; UDP galactose-4'-epimerase
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP200383) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)