

Product datasheet for **HK200173**

Glucokinase (GCK) Human qPCR Template Standard (BC001890)

Product data:

Product Type:	qPCR Template Standards
ACCN:	BC001890
Synonyms:	ATP:D-hexose 6-phosphotransferase; GK; GLK; glucokinase; glucokinase (hexokinase 4); hexokinase D, pancreatic isozyme; HHF3; HK4; HKIV; HXKP; MODY2; OTTHUMP00000024521; OTTHUMP00000159308
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP226189) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)