

Product datasheet for **HK200177**

GIMAP5 Human qPCR Template Standard (BC011732)

Product data:

Product Type:	qPCR Template Standards
ACCN:	BC011732
Synonyms:	FLJ11296; GTPase, IMAP family member 5; hIAN5; HIMAP3; IAN4; IAN4L1; IAN5; IMAP3; immune associated nucleotide 4 like 1; immunity associated protein 3
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP213051) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)