

## Product datasheet for **HK200178**

### GJB5 Human qPCR Template Standard (BC004379)

#### Product data:

Product Type:	qPCR Template Standards
ACCN:	BC004379
Synonyms:	connexin 31.1; CX31.1; gap junction protein, beta 5 (connexin 31.1); gap junction protein, beta 5, 31.1kDa; OTTHUMP00000004186
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP208437) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 <sup>7</sup> copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)