

## Product datasheet for **HK200185**

### GPR62 Human qPCR Template Standard (BC026357)

#### Product data:

Product Type:	qPCR Template Standards
ACCN:	BC026357
Synonyms:	G-protein coupled receptor; GPCR8; G protein-coupled receptor 62; KPG_005; MGC26943
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP216777) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 <sup>7</sup> copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)