

Product datasheet for **HK200188**

GPSM3 Human qPCR Template Standard (BC018724)

Product data:

Product Type:	qPCR Template Standards
ACCN:	BC018724
Synonyms:	C6orf9; G-protein signaling modulator 3 (AGS3-like, C. elegans); G-protein signalling modulator 3 (AGS3-like, C. elegans); G18; G18.1a; G18.1b; G18.2; NG1; OTTHUMP00000029352; OTTHUMP00000062706
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP214361) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)