

## Product datasheet for **HK200194**

### **HIST1H4J Human qPCR Template Standard (BC017361)**

#### **Product data:**

<b>Product Type:</b>	qPCR Template Standards
<b>ACCN:</b>	BC017361
<b>Synonyms:</b>	dj160A22.2; H4 histone family, member E; H4/e; H4F2iv; H4FE; histone 1, H4j; histone cluster 1, H4j; MGC29783; MGC166960
<b>Components:</b>	1 vial of lyophilized SybGREEN qPCR primer (HP214286) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 <sup>7</sup> copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
<b>Quality Control:</b>	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
<b>Storage:</b>	Store at -20°C.
<b>Stability:</b>	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)