

Product datasheet for **HK200202**

HTR1D Human qPCR Template Standard (BC007720)

Product data:

| | |
|------------------|---|
| Product Type: | qPCR Template Standards |
| ACCN: | BC007720 |
| Synonyms: | 5-HT1D; 5-hydroxytryptamine (serotonin) receptor 1D; HT1DA; HTR1DA; HTRL; OTTHUMP00000003270; RDC4 |
| Components: | 1 vial of lyophilized SybGREEN qPCR primer (HP200804) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer |
| Quality Control: | The kit has been tested to generate satisfactory qPCR data on ABI 7900HT |
| Storage: | Store at -20°C. |
| Stability: | The primer mix and template standard are stable for one year from date of shipping. |



[View online »](#)