

Product datasheet for **HK200215**

IL1 alpha (IL1A) Human qPCR Template Standard (BC013142)

Product data:

| | |
|------------------|---|
| Product Type: | qPCR Template Standards |
| ACCN: | BC013142 |
| Synonyms: | hematopoietin-1; IL-1A; IL1; IL1-ALPHA; IL1A (IL1F1); IL1F1; interleukin 1, alpha; preinterleukin 1 alpha; pro-interleukin-1-alpha |
| Components: | 1 vial of lyophilized SybGREEN qPCR primer (HP200543) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer |
| Quality Control: | The kit has been tested to generate satisfactory qPCR data on ABI 7900HT |
| Storage: | Store at -20°C. |
| Stability: | The primer mix and template standard are stable for one year from date of shipping. |



[View online »](#)