

Product datasheet for **HK200216**

IL24 Human qPCR Template Standard (BC009681)

Product data:

Product Type:	qPCR Template Standards
ACCN:	BC009681
Synonyms:	C49A; FISP; IL-24; IL-24 splice variant delE3; IL-24 splice variant delE5; IL10B; Interleukin 24; mda-7; MDA7; melanoma differentiation association protein 7; Mob-5; OTTHUMP00000034620; ST16; suppression of tumorigenicity 16 (melanoma differentiation)
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP209719) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)