

## Product datasheet for **HK200222**

### ITPKB Human qPCR Template Standard (BC015009)

#### Product data:

Product Type:	qPCR Template Standards
ACCN:	BC015009
Synonyms:	1D-myo-inositol-trisphosphate 3-kinase B; inositol 1,4,5-trisphosphate 3-kinase B; IP3 3-kinase; IP3K; IP3K-B; IP3KB; OTTHUMP00000035705; PIG37; proliferation-inducing protein 37
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP205955) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 <sup>7</sup> copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)