

Product datasheet for **HK200225**

KCNH6 Human qPCR Template Standard (BC006334)

Product data:

Product Type:	qPCR Template Standards
ACCN:	BC006334
Synonyms:	eag-related gene member 2; eag related protein 2; ERG2; ether-a-go-go related gene potassium channel 2; HERG2; Kv11.2; potassium voltage-gated channel, subfamily H (eag-related), member 6; potassium voltage-gated channel, subfamily H, member 6
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP215405) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)