

## Product datasheet for **HK200228**

### **BAT2L Human qPCR Template Standard (BC002872)**

#### **Product data:**

<b>Product Type:</b>	qPCR Template Standards
<b>ACCN:</b>	BC002872
<b>Synonyms:</b>	DKFZp781F05101; DKFZp781K12107; HLA-B associated transcript 2-like; KIAA0515; LQFBS-1; MGC10526; RP11-334J6.1
<b>Components:</b>	1 vial of lyophilized SybGREEN qPCR primer (HP210472) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 <sup>7</sup> copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
<b>Quality Control:</b>	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
<b>Storage:</b>	Store at -20°C.
<b>Stability:</b>	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)