

Product datasheet for **HK200234**

L2HGDH Human qPCR Template Standard (BC006117)

Product data:

Product Type:	qPCR Template Standards
ACCN:	BC006117
Synonyms:	2-hydroxyglutarate dehydrogenase; alpha-hydroxyglutarate oxidoreductase; alpha-ketoglutarate reductase; C14orf160; DURANIN; FLJ12618; L-2-hydroxyglutarate dehydrogenase; L-alpha-hydroxyglutarate dehydrogenase
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP215103) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)