

Product datasheet for **HK200260**

THRAP5 (MED16) Human qPCR Template Standard (BC017282)

Product data:

Product Type:	qPCR Template Standards
ACCN:	BC017282
Synonyms:	DRIP92; mediator complex subunit 16; THRAP5; thyroid hormone receptor-associated protein, 95-kD subunit; thyroid hormone receptor-associated protein 5; TRAP95
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP208607) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)