

Product datasheet for **HK200261**

MED18 Human qPCR Template Standard (BC002694)

Product data:

Product Type:	qPCR Template Standards
ACCN:	BC002694
Synonyms:	FLJ20045; mediator complex subunit 18; mediator of RNA polymerase II transcription, subunit 18 homolog; OTTHUMP0000003645; p28b; TRAP/mediator complex subunit p28b
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP234094) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)