

Product datasheet for **HK200263**

MFSD2 Human qPCR Template Standard (BC011587)

Product data:

Product Type:	qPCR Template Standards
ACCN:	BC011587
Synonyms:	FLJ14490; FLJ35904; FP1147; HMFN0656; major facilitator superfamily domain containing 2; OTTHUMP00000010893; PP10484; RP3-342P20.1
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP228651) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)