

Product datasheet for **HK200280**

Myeloid zinc finger 1 (MZF1) Human qPCR Template Standard (BC007777)

Product data:

Product Type:	qPCR Template Standards
ACCN:	BC007777
Synonyms:	myeloid zinc finger 1; MZF-1; MZF1B; Zfp98; zinc finger protein 42; zinc finger protein 42 (myeloid-specific retinoic acid-responsive); ZNF42; ZSCAN6
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP206943) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)