

Product datasheet for **HK200283**

NDUFB10 Human qPCR Template Standard (BC007509)

Product data:

Product Type:	qPCR Template Standards
ACCN:	BC007509
Synonyms:	NADH dehydrogenase (ubiquinone) 1 beta subcomplex, 10 (22kD, PDSW); NADH dehydrogenase (ubiquinone) 1 beta subcomplex, 10, 22kDa; NADH ubiquinone oxidoreductase PDSW subunit (RH 16p13.3); OTTHUMP00000158897; PDSW
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP207843) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)