

Product datasheet for **HK200292**

Neuropilin 1 (NRP1) Human qPCR Template Standard (BC007533)

Product data:

Product Type:	qPCR Template Standards
ACCN:	BC007533
Synonyms:	BDCA4; CD304; DKFZp686A03134; DKFZp781F1414; neuropilin-1; neuropilin 1; NP1; NRP; OTTHUMP00000020818; OTTHUMP00000020820; OTTHUMP00000020821; transmembrane receptor; VEGF165R
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP232738) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)