

Product datasheet for **HK200294**

HICE1 Human qPCR Template Standard (BC010176)

Product data:

Product Type:	qPCR Template Standards
ACCN:	BC010176
Synonyms:	Hec1-interacting and centrosome-associated 1; HEC1/NDC80 interacting, centrosome associated 1; MGC20533; NY-SAR-48; sarcoma antigen NY-SAR-48
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP227341) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)