

Product datasheet for **HK200303**

PDZD11 Human qPCR Template Standard (BC012996)

Product data:

Product Type:	qPCR Template Standards
ACCN:	BC012996
Synonyms:	AIPP1; ATPase-interacting PDZ protein; OTTHUMP00000023469; OTTHUMP00000023470; PDZ domain containing 11; PDZK11; PISP
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP212216) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)