

## Product datasheet for **HK200305**

### PIAS2 Human qPCR Template Standard (BC015190)

#### Product data:

|                  |   |
|------------------|---|
| Product Type:    | qPCR Template Standards   |
| ACCN:            | BC015190  |
| Synonyms:        | MGC102682; miz; MIZ1; Msx-interacting-zinc finger; PIASX; PIASX-ALPHA; PIASX-BETA; protein inhibitor of activated STAT, 2; protein inhibitor of activated STAT X; SIZ2; zinc finger, MIZ-type containing 4; ZMIZ4   |
| Components:      | 1 vial of lyophilized SybGREEN qPCR primer (HP207948) mix (1 nmol each primer, sufficient for 200 reactions)<br>1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 <sup>7</sup> copies, double-stranded DNA<br>1 vial containing 1 mL 10x PCR template standard dilution buffer |
| Quality Control: | The kit has been tested to generate satisfactory qPCR data on ABI 7900HT  |
| Storage:         | Store at -20°C.   |
| Stability:       | The primer mix and template standard are stable for one year from date of shipping.   |



[View online »](#)