

## Product datasheet for **HK200308**

### PLAC8 Human qPCR Template Standard (BC012205)

#### Product data:

Product Type:	qPCR Template Standards
ACCN:	BC012205
Synonyms:	C15; D5Wsu111e; onzin; OTTMUSP00000032100; OTTMUSP00000032101; OTTMUSP00000032102; placenta-specific 8
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP212306) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 <sup>7</sup> copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)