

## Product datasheet for **HK200310**

### PLEKHM2 Human qPCR Template Standard (BC008002)

#### Product data:

Product Type:	qPCR Template Standards
ACCN:	BC008002
Synonyms:	KIAA0842; novel RUN and PH domain-containing protein; pleckstrin homology domain containing, family M (with RUN domain) member 2; RP11-169K16.1; SifA (Salmonella-induced filaments A) and kinesin-interacting protein; SKIP
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP211437) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: $50 \times 10^7$ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)