

Product datasheet for **HK200312**

Adipose Triglyceride Lipase (PNPLA2) Human qPCR Template Standard (BC017280)

Product data:

Product Type:	qPCR Template Standards
ACCN:	BC017280
Synonyms:	1110001C14Rik; adipose triglyceride lipase; ATGL; DESNUTRIN; DKFZp667M109; FP17548; patatin-like phospholipase domain containing 2; PEDF-R; pigment epithelium-derived factor; transport-secretion protein 2.2; triglyceride hydrolase; TTS-2.2; TTS2
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP213638) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)