

Product datasheet for **HK200316**

PPP2R1A Human qPCR Template Standard (BC001537)

Product data:

Product Type:	qPCR Template Standards
ACCN:	BC001537
Synonyms:	alpha isoform of regulatory subunit A, protein phosphatase 2; medium tumor antigen-associated 61 KDA protein; MGC786; PP2A, subunit A, PR65-alpha isoform; PP2A, subunit A, R1-alpha isoform; PR65A; protein phosphatase 2 (formerly 2A), regulatory subunit A
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP210711) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)