

Product datasheet for **HK200330**

Rad51L1 (RAD51B) Human qPCR Template Standard (BC030219)

Product data:

Product Type:	qPCR Template Standards
ACCN:	BC030219
Synonyms:	DNA repair protein RAD51 homolog 2; hREC2; MGC34245; R51H2; RAD51 homolog B; RAD51-like 1; RAD51-like 1 (<i>S. cerevisiae</i>); RAD51B; REC2; RecA-like protein; recombination repair protein
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP226583) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)