

## Product datasheet for **HK200332**

### Retinoic Acid Receptor alpha (RARA) Human qPCR Template Standard (BC008727)

#### Product data:

Product Type:	qPCR Template Standards
ACCN:	BC008727
Synonyms:	NR1B1; nucleophosmin-retinoic acid receptor alpha fusion protein NPM-RAR long form; OTTHUMP00000164454; OTTHUMP00000164456; RAR; retinoic acid receptor, alpha; Retinoic acid receptor, alpha polypeptide
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP200901) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 <sup>7</sup> copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)