

Product datasheet for **HK200333**

RASGRF1 Human qPCR Template Standard (BC040275)

Product data:

Product Type:	qPCR Template Standards
ACCN:	BC040275
Synonyms:	CDC25; CDC25L; GNRP; GRF1; GRF55; guanine nucleotide-releasing factor, 55 kD; guanine nucleotide exchange factor; H-GRF55; PP13187; Ras-specific guanine nucleotide-releasing factor, CDC25 homolog; Ras-specific nucleotide exchange factor CDC25; Ras protein
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP206504) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)