

Product datasheet for **HK200348**

ROR gamma (RORC) Human qPCR Template Standard (BC031554)

Product data:

Product Type:	qPCR Template Standards
ACCN:	BC031554
Synonyms:	MGC129539; NR1F3; nuclear receptor ROR-gamma; OTTHUMP00000015253; RAR-re; RAR-related orphan receptor C, isoform a; RAR-related orphan receptor gamma; retinoic acid-binding receptor gamma; retinoid-related orphan receptor gamma; RORG; RZR-GAMMA; RZRG; TOR
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP230849) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)