

## Product datasheet for **HK200355**

### Serum Amyloid A (SAA1) Human qPCR Template Standard (BC007022)

#### Product data:

Product Type:	qPCR Template Standards
ACCN:	BC007022
Synonyms:	MGC111216; PIG4; SAA; serum amyloid A1; TP53I4; tumor protein p53 inducible protein 4
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP200312) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 <sup>7</sup> copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)