

Product datasheet for **HK200361**

SFRS17A (AKAP17A) Human qPCR Template Standard (BC028151)

Product data:

Product Type:	qPCR Template Standards
ACCN:	BC028151
Synonyms:	721P; B-lymphocyte surface antigen; CCDC133; CXorf3; DNA segment on chromosome X and Y (unique) 155 expressed sequence; DXYS155E; MGC39904; MGC125365; MGC125366; OTTHUMP0000022833; pseudoautosomal gene XE7; splicing factor, arginine/serine-ri; XE7; XE7Y
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP208291) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)