

Product datasheet for **HK200365**

SH2B1 Human qPCR Template Standard (BC010704)

Product data:

Product Type:	qPCR Template Standards
ACCN:	BC010704
Synonyms:	DKFZp547G11110; FLJ30542; KIAA1299; SH2 domain-containing putative adapter SH2-B; SH2-B; SH2-B alpha signaling protein; SH2-B gamma signaling protein; SH2-B homolog; SH2B; SH2B adaptor protein 1
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP234090) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)