

## Product datasheet for **HK200366**

### NSP3 (SH2D3C) Human qPCR Template Standard (BC032365)

#### Product data:

Product Type:	qPCR Template Standards
ACCN:	BC032365
Synonyms:	CHAT; FLJ39664; novel SH2-containing protein 3; NSP3; OTTHUMP00000022206; OTTHUMP00000022207; PRO34088; SH2 domain-containing 3C; SH2 domain containing 3C
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP230308) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 <sup>7</sup> copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)