

## Product datasheet for **HK200367**

### SLC19A1 Human qPCR Template Standard (BC003068)

#### Product data:

Product Type:	qPCR Template Standards
ACCN:	BC003068
Synonyms:	CHMD; FOLT; IFC1; reduced folate carrier; REFC; RFC1; solute carrier family 19 (folate transporter), member 1; solute carrier family 19 member 1
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP231560) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 <sup>7</sup> copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)