

## Product datasheet for **HK200368**

### SLC22A18AS Human qPCR Template Standard (BC030237)

#### Product data:

Product Type:	qPCR Template Standards
ACCN:	BC030237
Synonyms:	Beckwith-Wiedemann region 1B; Beckwith-Wiedemann syndrome chromosome region 1, candidate b; BWR1B; BWSCR1B; ORCTL2S; organic cation transporter-like 2 antisense; p27-BWR1B; SLC22A1LS; solute carrier family 22 (organic cation transporter), member 1-like an
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP209897) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 <sup>7</sup> copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)