

Product datasheet for **HK200380**

SOX10 Human qPCR Template Standard (BC002824)

Product data:

Product Type:	qPCR Template Standards
ACCN:	BC002824
Synonyms:	DOM; dominant megacolon, mouse, human homolog of; MGC15649; OTTHUMP00000028515; OTTHUMP00000195094; SRY (sex determining region Y)-box 10; SRY-related HMG-box gene 10; WS2E; WS4
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP209782) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)