

Product datasheet for **HK200383**

SPAG8 Human qPCR Template Standard (BC019247)

Product data:

Product Type:	qPCR Template Standards
ACCN:	BC019247
Synonyms:	BS-84; HSD-1; hSMP-1; MGC26201; OTTHUMP00000021352; SMP1; SPAG3; sperm associated antigen 8; sperm membrane protein 1
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP223555) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)