

Product datasheet for **HK200385**

Osteopontin (SPP1) Human qPCR Template Standard (BC017387)

Product data:

Product Type:	qPCR Template Standards
ACCN:	BC017387
Synonyms:	BNSP; bone sialoprotein I; BSPI; early T-lymphocyte activation 1; ETA-1; MGC110940; OPN; osteopontin; osteopontin/immunoglobulin alpha 1 heavy chain constant region fusion protein; secreted phosphoprotein 1 (osteopontin, bone sialopro; SPP1/CALPHA1 fusion
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP200549) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)