

Product datasheet for **HK200389**

BRDG 1 (STAP1) Human qPCR Template Standard (BC014958)

Product data:

Product Type:	qPCR Template Standards
ACCN:	BC014958
Synonyms:	BCR downstream signaling 1; BRDG1; OTTHUMP00000159640; signal-transducing adaptor protein-1; signal transducing adaptor family member 1; STAP-1
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP210109) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)