

Product datasheet for **HK200398**

TBRG4 Human qPCR Template Standard (BC014918)

Product data:

Product Type:	qPCR Template Standards
ACCN:	BC014918
Synonyms:	cell cycle progression 2 protein; CPR2; FASTKD4; FAST kinase domains 4; H_TD2522F11.8; KIAA0948; OTTHUMP00000159556; transforming growth factor beta regulator 4
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP208017) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)