

Product datasheet for **HK200400**

TCEB2 Human qPCR Template Standard (BC013306)

Product data:

Product Type:	qPCR Template Standards
ACCN:	BC013306
Synonyms:	ELOB; elongin, 18-kD subunit; elongin B; RNA polymerase II transcription factor SIII p18 subunit; SIII; SIII p18; transcription elongation factor B (SIII), polypeptide 2 (18kD, elongin B); transcription elongation factor B (SIII), polypeptide 2 (18kDa, el
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP226259) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)